Distributed by:

**Medtronic**

**Medtronic MiniMed**

18000 Devonshire Street  
Northridge, CA 91325  
USA  
800 646 4633  
818 576 5555

Medtronic B.V.  
Earl Bakkenstraat 10  
6422 PJ Heerlen  
The Netherlands

**Unomedical**  
A Convatec Company

**Legal Manufacturer:**  
Unomedical a/s  
Aaholmvej 1-3, Østed,  
4320 Lejre,  
Denmark  
Made in Mexico

Rx Only

€ € 0459

STERILE EO

REF MMT-100  
MMT-101


(01)00643169603523(10)NA

6025931-22A1\_b












## i-Port Advance™

Injection port • Port d'injection • Injektionsport  
Aplicador de inyecciones • Injectiepoort  
Porta di accesso per iniezioni • Injeksjonsport  
Injektioportti • Injektionsport • Injektionsport  
Porta de injeção • Injekcijski nastavek  
Injekční port • Injekciós bemenet • Port iniekcyjny  
Enjeksiyon portu • Injekčný port • Θύρα έγχυσης  
Инъекционный порт • 注射端口  
منفذ الحقن • فتحة التزجج

EN	Insertion Sites	CS	Místa zavedení
FR	Sites d'insertion	HU	Bevezetési helyek
DE	Einführstellen	PL	Miejsca wprowadzenia
ES	Zonas de inserción	TR	Giriş Bölgeleri
NL	Inbrengplaatsen	SK	Miesta zavedenia
IT	Siti di inserzione	EL	Σημεία εισαγωγής
NO	Innstikksteder	RU	Места введения
FI	Asetuskohdat	ZH	植入部位
SV	Appliceringsställen	AR	موقع الإدخال
DA	Indstikssteder	HE	מיקומות הזרה
PT	Locais de inserção	SL	Mesta vstavitev

EN	English	2
FR	Français	8
DE	Deutsch	14
ES	Español	20
NL	Nederlands	26
IT	Italiano	32
NO	Norsk	38
FI	Suomi	44
SV	Svenska	50
DA	Dansk	56
PT	Português	62
SL	Slovenščina	68
CS	Česky	74
HU	Magyar	80
PL	Polski	86
TR	Türkçe	92
SK	Slovenčina	98
EL	Ελληνικά	104
RU	Русский	110
ZH	中文	116
AR	128	العربية
HE	122	עברית

## **Indications for Use**

The i-Port Advance™ injection port is indicated for patients who administer or receive multiple daily subcutaneous injections of physician prescribed medications, including insulin. The device may remain in place for up to 72 hours to accommodate multiple injections without the discomfort of additional needle sticks. The i-Port Advance may be used on a wide range of patients, including adults and children.

## **Product Description**

The i-Port Advance is a subcutaneous injection port, which combines the injection port with an aid for insertion. The injection port is an integrated part of the inserter. Two injection port models are available: a 6 mm or 9 mm cannula length. The i-Port Advance is a sterile unit, delivered assembled for one time single use.

## **Contraindications**

- **Do not** use the same i-Port Advance for longer than 72 hours.
- **Do not** attempt to reapply the i-Port Advance after use.
- The i-Port Advance is neither intended nor indicated for use with blood or blood products.
- **Do not** use the i-Port Advance for any purpose other than the indications specified by your healthcare provider.
- **Do not** use a needle longer than 5/16 in. (8 mm) or thicker than 28 gauge when injecting into the i-Port Advance.

## **Warnings and Precautions**

- The i-Port Advance is only sterile and non-pyrogenic if the insertion device is unopened or undamaged. **Do not** use if the package is already open or is damaged, or if the protective liner has fallen off. Ensure the sterility by checking the sterile paper and tamper-proof seal.
- Carefully read through instructions prior to using the i-Port Advance. Failure to follow instructions may result in pain or injury.

- If the i-Port Advance is not securely placed in the inserter with the needle pointing straight ahead prior to insertion, pain or injury may occur.
- Be sure that the needle guard is removed before insertion.
- **Do not** use if introducer needle is bent or broken.
- **Do not** use a needle longer than 5/16 in. (8 mm). **Do not** attempt to use a 1/2 in. (12 mm) injection needle with the i-Port Advance. Doing so may damage the device causing unnecessary punctures of the skin, or tearing or puncturing of the soft cannula which may result in unpredictable medication delivery.
- **Do not** use a needle shorter than 3/16 in. (5 mm) when injecting into the i-Port Advance.
- **Do not** inject multiple drugs contraindicated for subcutaneous administration into the same site.
- The i-Port Advance is a single use device. Reuse of the port may damage the cannula in the device and increases the risk of the tape becoming loose. Reuse of the port may result in infection or site irritation and inaccurate medication delivery.
- Incorrect medication absorption, infection, or site irritation may result from improper application to the body or maintenance of the i-Port Advance insertion site, or both. If the insertion site becomes irritated or inflamed, remove and discard the device and apply a new device to a different location on the body.
- The soft cannula may crimp or bend when applying or while wearing. A crimped or bent cannula may limit or prevent medication delivery through the i-Port Advance. You should remove and discard the device immediately if you suspect a crimped or bent cannula. Indications that your i-Port Advance has a crimped or bent cannula include:
  - difficulty pressing down on the plunger of the syringe when injecting
  - injected medication escaping from the device onto the top of the septum
  - injected medication leaking onto the skin or saturating the tape
  - moisture present inside the transparent body of the device

- **Do not** reinsert the introducer needle into the i-Port Advance. Reinsertion could cause puncturing or tearing of the soft cannula which would result in unpredictable medication delivery or injury.
- Replace the i-Port Advance if the tape becomes loose or is displaced from its original site. Since the cannula is soft, it may not cause pain if it slips out, and this may occur without you noticing it. The soft cannula must always be completely inserted to receive the full amount of medication from the injection.
- Wash your hands and swab skin with alcohol or similar disinfectant where you intend to apply the i-Port Advance. Let skin surface dry completely prior to applying the device to the disinfected area.
- Ensure the septum is clean before every use (use pure water).
- **Do not** put perfumes or deodorants on the i-Port Advance as these may affect the integrity of the device.
- **Do not** inject into the i-Port Advance through clothing.
- Never point the loaded inserter towards any body part where insertion is not desired.
- Store i-Port Advance at room temperature. **Do not** store or leave the device in direct sunlight.
- Avoid applying the i-Port Advance to sites that contain scar tissue, lipodystrophy, or stretch marks. Avoid sites that are constrained by clothing or accessories, or that are subjected to rigorous movement during exercise. Avoid sites that are under a belt or on the waistline.
- **Do not** inject more than 75 times through a single device.
- The device has a residual volume up to 2.60 µL (0.0026 mL or 0.26 units). It is not known if the residual drug remains in the cannula or is absorbed into the subcutaneous tissue over time. The prescribing healthcare provider should consider if the retention of, or delay in delivery of small doses of a drug (such as 1 international unit or 10 µL of U100 insulin) may adversely affect


the therapeutic goals in some patients such as infants or small children.

## Recommendations

- Refer to your healthcare provider for instructions about insertion site rotation.
- Test your blood glucose per your healthcare provider's instructions to ensure that insulin delivery at the insertion site is appropriate.
- Test blood glucose 1 to 3 hours after inserting and using the i-Port Advance to ensure insulin delivery performs as anticipated, and measure your blood glucose on a regular basis.
- Avoid applying and using a new i-Port Advance prior to bedtime unless you can check your blood glucose 1 to 3 hours after injection.
- In the case of an unexplained high blood glucose, remove and discard the i-Port Advance and apply a new device to a different location on the body since the soft cannula could be dislodged, crimped or partially clogged. Should suggested actions not provide a resolution to the problem or should your blood glucose remain high, contact your physician or healthcare provider.
- After insertion, replace the lid onto the inserter to safely dispose of the inserter. For available sharps containers, please consult your local pharmacy.

# Instructions for Use

## i-Port Advance™

- |                        |                        |
|------------------------|------------------------|
| (A) Sterile liner      | (H) Lid                |
| (B) Center of inserter | (I) Round indentations |
| (C) Inserter           | (J) Septum             |
| (D) Finger grips       | (K) Tape               |
| (E) Introducer needle  | (L) Protective liner   |
| (F) Injection port     | (M) Soft cannula       |
| (G) Needle guard       |                        |

## Insertion

**Unfold the front cover of this booklet to view images that are numbered to the steps below.**

- 1** Wash your hands with soap and water.  
Select an insertion site on a recommended area (shown in gray).
- 2** Clean the insertion site with a disinfectant as directed by your healthcare professional. Be sure the area is dry before inserting i-Port Advance. If needed, remove hair around the insertion site to prevent improper attachment of the tape to the skin.
- 3** Pull the red tab and remove the seal.
- 4** Remove the sterile liner.

- 
- 5** Securely hold the inserter and remove the lid.
  - 6** Remove the protective liner from the tape.
  - 7** Twist and remove needle guard from introducer needle.
  - 8a** Place your fingers on the finger grips and hold the inserter so that the introducer needle is facing down.
  - 8b** With your other hand, pull up on the center of the inserter until it is locked in upright position.
  - 9** Place the inserter against your skin and place your fingers on the round indentations. Squeeze both round indentations at the same time to insert the injection port.
  - 10** Press down on the center of the inserter to secure the injection port to your skin.
  - 11** Grip the center of the inserter and gently pull to remove the inserter from the injection port.
  - 12** Press the tape securely to your skin.
  - 13** Put the lid back and dispose of the inserter according to local requirements.
- 

## Using i-Port Advance

- 
- 1** Clean the septum with pure water before each use.
  - 2** Puncture the septum using an injection needle.
  - 3** Slowly inject medication into the injection port.
  - 4** Remove the injection needle.
- 

## Removal

- 
- ⚠** Grip the edge of the tape and pull upward to remove the injection port from your skin.
- 

© 2015 Medtronic MiniMed, Inc. All rights reserved.

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Medtronic AB
Medtronic-Taiwan Ltd.
Medtronic (Thailand) Ltd.
Ticaret Ltd. Sirketi.
Medtronic B.V. Representative office in Ukraine
Medtronic Ltd.