


ENLITE™

RELEASED

Serter
User Guide

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
- F. Adhesive tab
- G. Sensor connectors
- H. Enlite overtape


1. Wash your hands.

WARNING: Wear gloves when inserting the sensor into someone other than yourself to avoid contact with patient blood. Minimal bleeding may occur.


2. Select an insertion site.


3. Clean the insertion site with alcohol. Let the area air dry.


4. Open the sensor package.


5. Remove the pedestal from the package. Place the pedestal on a flat surface.


6-7. Push the Serter down onto the pedestal until the base of the Serter sits flat on the table.

8. To detach the sensor from the pedestal, place two fingers on the pedestal arms and slowly pull the Serter straight up. Do not detach the pedestal from the Serter in mid-air, as this might damage the sensor.

WARNING: Never point a loaded Serter toward any body part where insertion is not desired. An accidental button-push may cause the needle to inject the sensor in an undesired location causing minor injury.


9. Place the base of the Serter flat against your insertion site.

WARNING: The sensor insertion process is a multi-step process (as shown below). Not following directions may result in an accidental needle stick, incorrect sensor insertion, or the sensor being pulled out of the body. When following the steps below, do not remove the Serter until you have pressed (10a), released (10b), and then pressed the button again (10c) to ensure correct sensor insertion. The Serter injects the sensor upon the first button press and release. The second button press is required to detach the Serter from the sensor.


Note: After you insert the sensor (10a-10c), hold the button in while lifting the Serter off the sensor insertion site (10d).

To insert the sensor, complete the following four steps:


10a

10a. Press the button in and then release.

10b

10b. Wait 5 seconds to allow the adhesive time to stick to skin.



10c

10c. Press and hold in the button again. Do not release it.

10d

10d. While holding the button in, slowly lift the Serter away from the skin.

If you inserted the sensor into yourself, complete step 11a. If you are a healthcare professional or caregiver who inserted the sensor into a patient, complete step 11b.

Patient:


11a. Gently hold the base of the sensor against the skin. Hold the needle housing at the top and slowly pull straight up, away from the sensor.

OR

Healthcare professional or caregiver:


11b. Wrap sterile gauze around the sensor (as shown in image 11b). Gently hold the base of the sensor against the skin. Hold the needle housing at the top and slowly pull straight up, away from the sensor.

WARNING: Watch for bleeding at the insertion site. If bleeding occurs at your sensor site (under, around, or on top of the sensor), apply steady pressure using sterile gauze or a clean cloth placed on top of the sensor for up to three minutes. The use of unsterile gauze can cause an infection. If bleeding does not stop, then remove the sensor and apply steady pressure until the bleeding stops.



12. Hold the sensor in place and gently remove the white paper from under the adhesive pad. Firmly press the adhesive against the skin for a minimum of five seconds.

NOTE: Medtronic adhesives are pressure sensitive. Apply pressure for a minimum of five seconds to ensure that the sensor remains adhered to the skin throughout the wear period.


13. Straighten the sensor adhesive tab so that it lays flat against the skin.



14. Remove the protective white paper from the overtape backing. Do not remove the white paper from the sides of the overtape.

15. Apply the adhesive to the rounded part of the sensor and the skin.



16. Stretch the overtape around the connectors so that it lays flat against the sensor adhesive pad.

17. Remove the protective white paper from the sides of the overtape and smooth the adhesive to the skin.

Example


These images show examples of the Enlite overtape applied correctly.

18

MiniMed® 530G System
User Guide


18. Refer to the *Sensor features* chapter, found in your insulin pump user guide, for details on how to enter your sensor settings into your pump.

19


19. Connect the transmitter to the sensor.

20



20. Wait for the transmitter to flash.

21-23. Cover the transmitter with the sensor adhesive tab.

CAUTION: Do not pull the adhesive tab on the sensor too tight across the transmitter. Doing so may prevent a proper connection between the transmitter and sensor, or may cause the sensor to pull out of your body.

RELEASED

RELEASED

Introduction

The Enlite™ insertion device (Serter) is used as an aid for the insertion of a Medtronic Diabetes glucose sensor.

Potential risks related to Serter use

General risks with Serter use may include skin infection around the area where the Serter is used.

Indications for use

The Enlite Serter is indicated for use for a single patient as an aid to insert the Enlite glucose sensor (sensor).

Contraindications

Do not use the Enlite Serter on products other than the Enlite sensor. Medtronic cannot guarantee this product's safety or efficacy if used on other products.

General warnings

Read this entire user guide before attempting to insert the sensor. The Enlite Serter does not work the same as other Medtronic insertion devices. **The Serter injects the sensor into the insertion site when the button is released, not when the button is pressed.** Failure to follow directions may result in pain or injury.

Do not make therapy treatment decisions based on sensor glucose values because sensor glucose and blood glucose values may differ. If your sensor glucose reading is low or high, or if you feel symptoms of low or high glucose, confirm your blood glucose with your BG meter prior to making therapy decisions to avoid severe low or high glucose conditions.

Keep an arm's length between you and the patient when removing the needle. A retractable needle is attached to the sensor and minimal blood splatter may occur.

Keep the needle housing within sight at all times to avoid accidental needle stick injury.

Always inspect the sensor packaging for damage prior to use. Sensors are sterile and nonpyrogenic, unless the package has been opened or damaged. Do not use any of the sensors if the sterile package has been opened or damaged. Use of an unsterile sensor can cause site infection.

Watch for bleeding at the insertion site (under, around, or on top of the sensor). Apply steady pressure using sterile gauze or a clean cloth placed on top of the sensor for up to three minutes, if bleeding occurs. The use of unsterile gauze can cause an infection.

Do not attach the transmitter to the sensor, if bleeding does not stop. Blood can get into the transmitter connector, which may damage the device. Remove the sensor and continue to apply steady pressure until the bleeding stops.

1. Remove the sensor and discard in a sharps container.
2. Check the site for redness, bleeding, irritation, pain, tenderness, or inflammation and treat accordingly.
3. Insert a new sensor in a different location.

Make sure the sensor is securely placed in the Serter to avoid pain or minor injury.

Keep the Serter away from children. This product contains small parts and may pose a choking hazard.

Refer to the Enlite Glucose Sensor User Guide for additional sensor related warnings.

General precautions

Wash your hands with soap and water before inserting the sensor to help prevent site infection.

Wear gloves when inserting the sensor into someone other than yourself to avoid contact with patient blood.

Do not insert the sensor through tape to ensure proper sensor function and insertion.

Only use alcohol to prepare the insertion site to ensure that residue is not left on the skin.

Rotate the sensor insertion site so that sites do not become overused.

Discard used sensors and needle housings in a sharps container after each use to avoid accidental needle stick injury.

Do not clean, resterilize, or try to extract the needle from the needle housing. An accidental needle stick may occur.

Do not reuse sensors. Reuse of the sensor may cause damage to the sensor surface and lead to inaccurate glucose values, site irritation, or infection.

Where to insert the sensor

CAUTION: Avoid the 2 inch (5.0 cm) around the navel to help ensure a comfortable insertion site, and to help with adhesion.

Choose an insertion site that has an adequate amount of subcutaneous fat. The Enlite sensor has been approved for use on the abdomen. See step 2 for a diagram of these areas.

Do not insert the sensor in muscle or areas constrained by clothing or accessories; areas with tough skin or scar tissue; sites subjected to rigorous movement during exercise; or in sites on the belt or waistline to avoid chafing or site irritation.

Storage and handling

Store the Serter in the released position, at room temperature, to maintain optimum product performance and life.

Specifications

Dimensions	2.93 x 1.95 x 2.45 inches	Weight	71.2 grams
------------	---------------------------	--------	------------

Serter life of use

The Serter can be used for up to 600 sensor insertion uses. And, it has been tested through 92 disinfection cycles for functionality.

WARNING: The disinfection process may damage the Serter prior to 600 uses. Always inspect the Serter for any signs of cracking, flaking, or damage prior to use. Discard the Serter if you see any of these signs to avoid incorrect sensor insertion, pain or injury.

Serter maintenance

The Serter is intended for single-patient use. Always clean and disinfect the Serter after each use.

To clean and disinfect the Serter, you will need to following materials: gloves, mild liquid soap (for example Ivory® liquid soap), household bleach (for example Clorox® Regular-Bleach), soft-bristled toddler toothbrush, and a container.

To clean the Serter:

1. Prepare a mild liquid soap solution using one (1) teaspoon of mild liquid soap per one (1) gallon of tap water.
2. Rinse the Serter under room temperature tap water for at least one minute. Continue rinsing until the Serter is visibly clean. Make sure all hard-to-reach areas are rinsed completely.
3. While rinsing the underside of the Serter, press and release the mechanism on the underside of the Serter to ensure that the underside of the Serter is completely rinsed with water.
4. Submerge the Serter in the mild liquid soap solution and soak for at least 10 minutes.
5. Brush the entire surface of the Serter using a soft-bristled toddler toothbrush until visibly clean. Pay close attention to brushing the hard-to-reach areas.
6. Rinse the Serter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.
7. Shake off any excess water or moisture, and place the Serter upright on a clean, dry, cloth to air-dry.

To disinfect the Serter, continue with steps 8 through 11:

8. While wearing gloves, prepare a 1:10 bleach solution by using one (1) cup of household bleach to nine (9) cups of room temperature tap water.
9. Submerge the Serter in the bleach solution and soak for 30 minutes.
10. Rinse the Serter under room temperature tap water for three minutes, ensuring all hard-to-reach areas are rinsed completely.
11. Shake off any excess water or moisture, and place the Serter upright on a clean, dry, cloth to air-dry.




Icon table

	Follow instructions for use.	REF	Catalogue number
	Batch code		Manufacturer

©2012, Medtronic MiniMed, Inc. All rights reserved. Enlite™ is a trademark of Medtronic MiniMed, Inc. MiniLink® is a registered trademark of Medtronic MiniMed, Inc. Clorox® is a registered trademark of The Clorox Company. Ivory® is a registered trademark of The Procter & Gamble Company.

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